<b>de</b>	Gebrauchsanweisung 📄 4	<b>pl</b>	Instrukcja obsługi 📄 48
<b>en</b>	Instructions for Use 📄 8	<b>ru</b>	Руководство по эксплуатации 📄 52
<b>fr</b>	Notice d'utilisation 📄 12	<b>hr</b>	Upute za uporabu 📄 56
<b>es</b>	Instrucciones de uso 📄 16	<b>sl</b>	Navodilo za uporabo 📄 60
<b>pt</b>	Instruções de utilização 📄 20	<b>sk</b>	Návod na použitie 📄 64
<b>it</b>	Istruzioni per l'uso 📄 24	<b>cs</b>	Návod na použití 📄 68
<b>nl</b>	Gebruiksaanwijzing 📄 28	<b>bg</b>	Ръководство за работ 📄 72
<b>da</b>	Brugsanvisning 📄 32	<b>ro</b>	Instrucțiuni de utilizare 📄 76
<b>fi</b>	Käyttöohjeet 📄 36	<b>hu</b>	Használati útmutató 📄 80
<b>no</b>	Bruksanvisning 📄 40	<b>el</b>	Οδηγίες Χρήσης 📄 84
<b>sv</b>	Bruksanvisning 📄 44	<b>tr</b>	Kullanma talimatları 📄 88

**Not for sale in the US market**

## Dräger DrugCheck® 3000 IVD

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## 1 For your safety

### 1.1 Follow the Instructions for Use

Any use of the Dräger DrugCheck® 3000<sup>1)</sup> requires complete understanding and strict observation of these Instructions for Use. The Dräger DrugCheck 3000 is only approved for professional application in in-vitro diagnostics.

The components of the Dräger DrugCheck 3000 are only intended for the described use.

## 2 Description

### 2.1 Product overview

#### Dräger DrugCheck 3000 Test Kit (fig. A)

- 1 Sample collector with colour indicator
- 2 Test cassette with indicator ring
- 3 Funnel-shaped test opening
- 4 Security seal
- 5 result window with test strip

#### Result window with test strip (fig. B)

- 6 Control line area (C)
- 7 Test line area (MET, AMP, THC, BZO, OPI, COC)

Each packing unit contains:

- 20 x Dräger DrugCheck 3000 test kits, each packed in a sealed foil pouch
- 1 x Instructions for Use

The Dräger DrugCheck 3000 test kit consists of 2 components: the test cassette and the sample collector

The removable sample collector has a red colour indicator to display sufficient sampling.

The test cassette contains an ampoule with buffer solution, a red coloured indicator ring for handling control, as well as 2 test strips with 2 control lines (C) and up to 6 strips for the various drugs/drug classes, depending on the product variant.

#### Additional recommended material

- Disposable gloves (e. g. latex or nitrile gloves)
- Timer

### 2.2 Intended use

The Dräger DrugCheck 3000 is a point of collection test for the simultaneous detection of up to 6 drugs/drug classes (depending on the product variant) in human oral fluid for diagnostic purposes (in-vitro diagnostic device) and/or forensic use. These drugs/drug classes are cocaine, opiate, amphetamine, methamphetamine, THC/cannabis and benzodiazepine. Each drug/drug class is represented by a separate test line on one of the test strips in the result window on the test cassette.






The Dräger DrugCheck 3000 packaging gives all important information on the drugs/drug classes with their abbreviations (see also the table on "2.2.1"), detection limits (in ng/mL), shelf life and storage conditions.

The detection limit (cut-off) is the deciding factor regarding the presence of the drugs/drug classes in the sample. The test result is correspondingly positive or negative above or below the limit value. To obtain a confirmed analytical result a more specific method must be used. The preferred method is the liquid or gas chromatography/mass spectrometry (GC-MS or LC-MS). All the results from the Dräger DrugCheck 3000 require professional assessment with reference to further clinical evaluation of the test subjects. This applies in particular to a positive ("not negative") test result, see also table 1 (T1) on page 92.

### 2.2.1 Drug (target analyte)

Drug	Target analyte
AMP	Amphetamine
BZO	Benzodiazepine
COC	Cocaine
MET	Methamphetamine
OPI	Opiates
THC	Tetrahydrocannabinol

### 2.3 Explanation of symbols

	Strictly follow the Instructions for Use!
	Manufacturer
	Single use only
	In-vitro diagnostic medical product
	Batch
	Temperature limit during storage
	Expiry date
	Contents sufficient for <n> tests
	Measurement mode: fast (THC)
	Measurement mode: sensitive (THC)

### 2.4 Test principle

The Dräger DrugCheck 3000 is based on the immunoassay principle of competitive inhibition. Drugs contained in the sample compete with drugs on the test membrane for the binding of nanoparticles coated with antibodies.

#### Sampling:

The Dräger DrugCheck 3000 is designed for use with saliva samples, which are taken with the integrated sampler.

The sample does not require special treatment. The sample is collected directly into the sample collector via absorption. The sample collector is then inserted into the funnel-shaped test opening in the test cassette to perform the test.

#### Test:

The fluid sample is then rinsed in the buffer solution by shaking. After waiting for the pre-incubation time to elapse, the test can be started by breaking the security seal and pressing the sample collector downwards with the test cassette.

Gold particles coated with antibodies react with drug molecules from the fluid sample and during the course of the test with drug conjugates on the test membrane. If the sample is drug-free, the antibodies can react freely with the drug conjugates, forming a red line on the test strip. If the drug is available in sufficient concentration, the gold particles are impeded from adhering to the drug conjugate. If the sample is preliminarily positive ('not negative'), no test line will appear in the area for the corresponding drug in the test strip result window.

#### Quality control:

An additional function check independent of the fluid sample is integrated into every test membrane. The appearance of a red line in area C in the viewing window indicates that sufficient saliva was collected and the test membrane wetted fully and that the antibody/antigen reaction occurred as expected.

1) Dräger DrugCheck is a registered trademark of Dräger.

## 2.5 Precautions and warnings

### NOTICE

For hygiene reasons, use gloves when handling the Dräger DrugCheck 3000. Do not touch the test cassette with bare hands before and after sampling and follow the current hygiene rules.

- Follow proper handling and disposal procedures.
- Open the foil pouch only immediately before use to prevent contamination of the sampler.
- Do not use the Dräger DrugCheck 3000 if the foil pouch is damaged in any way.
- Do not use the Dräger DrugCheck 3000 beyond the expiration date on the foil pouch. The expiry date is in the format YYYY-MM. Example: 2019-01 means that the Dräger DrugCheck 3000 is not permitted to be used after the end of January 2019.

## 3 Testing

### 3.1 Preparing the test

### NOTICE

Results may be incorrect if maximum storage times are not observed.

- Ensure that the test subject has not eaten or drunk anything for at least 10 minutes before sampling, e. g. chewing gum, tobacco, coffee or other substances.
- Ensure that the ambient temperature is between 5 °C and 40 °C and the relative air humidity is between 5 % and 95 %.
- Open the foil pouch.
- Only remove the sample collector immediately before using the test kit.

### 3.2 Test procedures

#### 3.2.1 Collecting the sample

1. Remove the sample collector from the test kit (**fig. 1**) and hand to the test subject.
2. Instruct the test subject to keep the sample collector in the cheek pouches and move it carefully from one side of the mouth to the other (**fig. 2**). This process takes approx. 15 seconds. If the red colour indicator on the sample collector mouthpiece loses its colour, the sampling process can be ended. The test subject must not suck or chew on the sample collector!
3. Insert the sample collector into the funnel-shaped test opening in the test cassette and press down as far as it will go to break the ampoule (**fig. 3**). There must not be any gap between the sample collector and test cassette.

#### 3.2.2 Preparing and starting the test

4. Shake the test kit **well** until the red indicator ring has lost all colour (**fig. 4**). This process takes at least 15 seconds. A small amount of residual colour in the indicator does not affect the evaluation.
5. Place the test kit on an even surface (**fig. 5**) or hold upright in your hand. Depending on the required detection limit (cut-off) for the detection of cannabis, the following waiting times (pre-incubation) must be observed:
  - a. Measurement mode "sensitive" (THC)  Wait **60 seconds**, then proceed.
  - b. Measurement mode "fast" (THC)  Wait **10 seconds**, then proceed.
6. After the waiting time has elapsed, break the security seal on the test kit and break off (**fig. 6.1**).
7. Push the test kit in an upright position at the sample collector forcefully down to the limit stop (**fig. 6.2**) and tap the base lightly on a flat surface (**fig. 6.3**). The saliva sample now flows out along the test strip.
8. If the saliva sample does not flow along the test strip within one minute, tap the base of the test kit on a flat surface again. Continuously monitor the lines forming on the test strip.

## 4 Evaluation

### 4.1 Requirements for the evaluation

Ensure that the lighting is adequate when reading off the results (e. g. daylight or torch).

### 4.2 Test results

#### Control lines

The control lines form in the top section C of the view window of the test strip after starting the test.

#### Negative test result

Once test lines appear in the drug/drug classes area in the result window and the control lines appear as well, the "negative" results can be read off for the corresponding drug/drug class. This means there are no drugs above the detection limit in the sample. The intensity of the lines can vary, meaning that weak, partially coloured or broken lines should also be evaluated as "negative".


#### Positive test result

If no red line appears in the drugs/drug classes area of the result window for **5 minutes** after starting the test, with a control line present, this indicates a positive test result for this drug/drug class. This means there are drugs above the detection limit in the sample.

#### Test valid/invalid

The test results must be evaluated within **10 minutes** after starting the test (for start time, see **fig. 6.2**). If the control lines do not appear on one or both test strips within 10 minutes after the start time of the test, the corresponding test strip must not be evaluated and is invalid. The intensity of the control lines may vary. For this reason, only the complete lack of a red line in area C of the result window for the corresponding test strip is a sign of an invalid test. Repeat the analysis if a statement is required on the invalid drug/drug class.

#### Examples for the evaluation (variant for 6 drugs/drug classes):

The test is valid and negative.

The test is valid and positive for opiates.

The test is only valid for the left test strip and it is negative. The right test strip is invalid.

### NOTICE

With immuno-chemical tests, factors may occur that affect the reliability. In rare cases this means that with a positive result the test subject has not always actually consumed drugs (false positive result). For this reason, a more comprehensive confirmation analysis (e.g. GC-MS or LC-MS) is recommended to confirm all positive test results. See chapter 6 "Limitations of the procedure".

## 5 Quality control

Regular quality control is part of good analytical practice and may be required by the relevant authorities. Always contact the responsible authorising and accreditation bodies to ensure that the quality control programme implemented corresponds to the applicable standards.

#### Internal control

An integrated process control is conducted with every analysis (see also chapter 2.4 "Test principle"). If the operation has been carried out without errors and the result can be reliably evaluated, the control line forms in area C of the result window. If this line does not appear, the test must be evaluated as invalid and not applicable. The test should then be repeated.

#### External control

The features of the DrugCheck 3000 allow quality control to be run on site by taking a drug-free sample and processing it. The result for all drugs/drug classes must be negative in this case. When using positive controls ensure that they are suitable for the Dräger DrugCheck 3000.

## 6 Limitations of the procedure

The analytical features can deviate at the limit of the operating temperature range. At temperatures below +10 °C, positive test results must not be evaluated until 10 minutes after starting the measurement. We recommend warming the test kit in your hands before use.

Start the test kit between -5 °C and 5 °C as described in Section 3.2.2 and move it to a warm environment (e.g. a warm car). The test kit can also be placed in its foil pouch and inserted upright in a jacket pocket. Grasp the foil pouch with the test kit inside with a warm hand to warm it up. Read the test results 10 minutes after the measurement was taken.

At temperatures over +35 °C and if there is a control line, positive test results must be evaluated 3 minutes after starting the test.

Due to non-specific interactions (physiological deviation, state of health or sample contamination), false positive or false negative results may occur in rare cases.

The antibodies used in the Dräger DrugCheck 3000 have been developed specifically for the detection of illegal substances (drugs/drug classes). For the user, this means for example that the antibody reagents of the Dräger DrugCheck 3000 are reacting with chemically similar substances such as prescription drugs or non-prescription medicines and are resulting in false positive results. The test has objectively detected the chemically related substance, in the confirming lab analysis, a corresponding illegal drug cannot be proven, however.

To achieve a confirmed analytical result, an alternative chemical method such as GC-MS or LC-MS must be used. All the results of the Dräger DrugCheck 3000 require professional assessment with reference to further clinical evaluation of the test subjects. This applies especially for positive results.

The sampling process can be closely monitored. This means that faking the sample is unlikely. However, if faking or substitution of a sample is suspected, dispose of the sample and repeat the test with a new Dräger DrugCheck 3000 test kit.

## 7 Performance features

### 7.1 Influence of food and beverages

Saliva may be contaminated immediately after consuming food, beverages or confectionery. To determine any influence of the contamination on the results of the Dräger DrugCheck 3000, saliva samples were taken and interpreted immediately after consuming sample preparations:

cola, toothpaste, American cranberry juice, aseptic mouthwash, water, cough medicine (not containing codeine), coffee, chewing gum, chocolate, fruit tea, herbal tea, cigarettes, orange juice, full fat milk, cough sweets containing menthol, mouthwash and tomatoes.

None of the interpretations yielded incorrect results.

As it is impossible to detect the potential influence of all foods on the test, it is strongly recommended to observe a waiting time of **10 minutes** before sampling.

## 7.2 Analytical features

### 7.2.1 Analytical specificity

Most of the immuno-chemical detection reactions are not monospecific for a single analyte (drug) but they react to a group of analytes (drug class) with a similar chemical structure (e.g. the opiate test of the Dräger DrugCheck 3000 detects a number of different opiates).

Therefore we do not recommend that you use the results of a Dräger DrugCheck 3000 as a basis for (semi) quantitative statements on the concentrations of one single analyte or a group of analytes. Data on the analytes which can be detected with the corresponding test of the DrugCheck 3000 as well as data on the individual concentrations which separately generate a positive result can be found in table 1 (T1) on page 92. The analytes found in table 2 (T2) on page 92 are not detected at concentrations below **10,000 ng/mL** by the Dräger DrugCheck 3000.

### 7.2.2 Repeatability

Reproducibility studies were conducted with commercially available reference standards and negative saliva samples. Every saliva sample was enriched with corresponding standards to maintain the desired concentration of the test analytes (no drugs, 200 % limit value). Every sample was tested at every analyte concentration ten times over three different days with the same batch of the Dräger DrugCheck 3000.

Correspondence with the expected test results was achieved in 100 % of the cases.

## 7.3 Analytical features

To test the analytical features of the Dräger DrugCheck 3000, native saliva samples were collected and analysed using the Dräger DrugCheck 3000. Some of the native saliva samples were enriched using forensic drug standards at defined concentration values. The analytical features of the Dräger DrugCheck 3000 are given in table 3 (T3) on page 92. The specified analytical performance is based on tests with the specified detection limits (cut-off).

## 8 Notes on troubleshooting

Possible fault	Possible fault cause	Preventative measures
No control line appears.	No or low colouration of the indicator ring. Not shaken for long enough or strongly enough.	Follow the Instructions for Use.
Sample does not completely flow up the test strip.	Dräger DrugCheck 3000 is not positioned on an even surface. Dräger DrugCheck 3000 has not been started when in a vertical position. The test strip has no contact with the sample liquid.	Follow the Instructions for Use. Observe the waiting time.  Lightly tap the bottom of the test kit repeatedly on a flat surface.
Unclear results or blurred lines	The test strips have been soaked with too much of the sample. Dräger DrugCheck 3000 has not been correctly started. Improper positioning of the Dräger DrugCheck 3000 Dräger DrugCheck 3000 is not positioned on an even surface.	Follow the Instructions for Use.
Colour blindness (for analyte result interpretation)	The result and control lines are coloured. The intensity of the colour is not significant as only the existence or absence of a line is relevant for the interpretation of the result.	Follow the Instructions for Use.
Questionable results	Physical degeneration of the Dräger DrugCheck 3000, incorrect storage, opening the packaging too early before starting the test or reading off the results too late.	Follow the Instructions for Use for the correct device storage, handling and interpretation of the results.
	Multiple occurrence of invalid results, sample is contaminated.	Observe waiting time of 10 minutes before starting the sampling process.
Formation of test and control lines is delayed.	Temperature too low.	Warm the Dräger DrugCheck 3000 in your hands before use.  Wait for the runtime of 10 minutes.
	Indicator ring insufficiently discoloured.	Shake the test kit until the indicator ring is clearly discoloured.
False-positive results at low temperatures.	The test lines are evaluated too early.	Wait for the runtime of 10 minutes.

## 9 Storage

The Dräger DrugCheck 3000 must be stored at temperatures between +4 °C and +25 °C.

## 10 Disposal

The Dräger DrugCheck 3000 test kits can be disposed of in domestic waste.

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## T1 – Specificity

Cocaine related compounds	[ng/mL]
Benzoylcegonine	400
Cocaethylene	500
Cocaine	20
Ecgoninemethylester	10000
Norcocaine	100
Procaine	10000

Opiates related compounds	[ng/mL]
6-Monoacetylmorphine	35
Buprenorphine	10000
Codeine	20
Dihydrocodeine	10
GHB	10000
Hydrocodone	10
Hydromorphone	10
Morphine	20
Morphine-3β-D-glucuronide	1000
Nalorphine	35
Naloxone	10000
Norcodeine	10000
Normorphine	100000
Oxycodone	1000

THC related compounds	[ng/mL]
Cannabidiol	90000
Cannabinol	350
delta8-THC	5000
THC (delta9-THC)	15 (25)
THC-COOH (11-Nor-9-Carboxy-delta9-THC)	50
THC-OH (11-Hydroxy-delta9-THC)	50
CP 47,797	>100000
JWH-18	>100000

Amphetamine related compounds	[ng/mL]
Dopamine	40000
MBDB	50000
MDA	50
MDEA	10000
Phentermine	100000
S-(+)-Amphetamine	35
Tyramine	5000

Methamphetamine related compounds	[ng/mL]
Chloroquine sulphate	5000
Ephedrine	1000
Fenfluramine	300
MBDB	70
MDA	10000
MDEA	10000
MDMA	75
Procaïn	4000
Pseudoephedrine	100000
R-(+)-Methcatinone	2000
Ranitidine	20000
S-(-)-Methcatinone	4500
S-(+)-Methamphetamine	35

Benzodiazepine related compounds	[ng/mL]
Aminoflunitrazepam	200
Alprazolam	10
Bromazepam	25
Chlordiazepoxide	3000
Clonazepam	60
Desalkylflurazepam	25
Diazepam	15
Flunitrazepam	20
Flurazepam	>10000
Lorazepam	100
Midazolam	80
Nitrazepam	10
Nordiazepam	15
Oxazepam	30
Prazepam	>10000
Temazepam	5
Triazolam	20

## T2 – No Interferents

Substance		
Acetaminophen	Acetylsalicylic acid	Ascorbic acid
Caffeine	Cotinine	Ibuprofen
Imipramine	Naproxen	Nicotine
Phenobarbital	PCP	Tetracycline
Tramadol		

## T3 – Analytical Performance

Cocain:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<10	10-20	>20-30	>30
Positive	0	0	10	40
Negative	140	10	0	0

Agreement 100 %

Opiates:  
Morphin:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<10	10-20	>20-30	>30
Positive	0	0	10	40
Negative	140	10	0	0

Agreement 100 %

Tetrahydrocannabinol:  
(-)-Delta-9-Tetrahydrocannabinol:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<7.5	7.5-15	>15-22.5	>22.5
Positive	2	0	7	40
Negative	138	10	3	0

Agreement 97.5 %

Amphetamine:  
S-(+)-Amphetamine:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<17.5	17.5-35	>35-52.5	>52.5
Positive	1	0	10	40
Negative	140	10	0	0

Agreement 100 %

---

Methamphetamine:  
S-(+)-Methamphetamine:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<17.5	17.5-35	>35-52.5	>52.5
Positive	0	2	10	40
Negative	140	8	0	0

Agreement 99 %

Benzodiazepine:  
Diazepam:

Dräger DrugCheck 3000	Reference Sample Contains [ng/mL]			
	<7.5	7.5-15	>15-22.5	>22.5
Positive	0	1	10	40
Negative	140	9	0	0

Agreement 99.5 %

 Directive 98/79/EC

**Dräger Safety AG & Co. KGaA**

Revalstraße 1  
23560 Lübeck, Germany  
Phone+49 451 882 0  
Fax +49 451 882 20 80  
[www.draeger.com](http://www.draeger.com)

**Dräger Safety AG & Co. KGaA**

Revalstraße 1  
23560 Lübeck, Germany  
Phone +49 451 882 0  
Fax +49 451 882 20 80  
[www.draeger.com](http://www.draeger.com)

**90 33 736** - GA 4755.210\_MUL077  
© Dräger Safety AG & Co. KGaA  
Edition 06 - April 2018  
Subject to alteration